

CD-3A CARBON DIOXIDE ANALYZER



The CD-3A Carbon Dioxide Analyzer provides a continuous and accurate measurement of gases containing up to 15% CO₂. The CD-3A has the rapid response necessary for breath-by-breath analysis during respiration. It is also ideal for laboratory applications requiring high accuracy and sensitivity.

Exclusive Features:

- Wide range of Gas Detection 0 to 15% CO₂
- Superior accuracy ± 0.02 carbon dioxide
- Very high sensitivity $\pm 0.001\%$ over the entire 0 - 15% CO₂ range
- Extremely rapid 25 millisecond response - needed for accurate breath-by-breath measurements
- End Tidal CO₂ Displayed - at the touch of a button
- Reliable, low maintenance operation

Exact VO₂, VCO₂ & RER Measurements are too important for speculation, choose the Gas Analyzers researchers have depended on for over 45 years.

CONTACT US at SALES@AEITECHNOLOGIES.COM or + 1 630 548 3545 for a **QUOTE**
visit us at www.aeitechnologies.com