




TECHNOLOGIES

Cardio Stress ECG

Cardio Card: PC based Stress ECG (Stress PC ECG) / PC EKG (Electrocardiograph) System
(Desk-top or NoteBook Formats)

- Simultaneous PC Based ECG (12 Lead PC ECG) / PC EKG Diagnostic Acquisition.
- WINDOWS: Point & Click User Interface.
- Real-Time 3, 6 or 12 Lead Color Display of ECG / EKG Complex.
- 5-15 ft. Patient ECG / EKG Cable w/ Locking Replaceable Leads (plus 3 ft.).
- Full Patient Demographics & Information.
- HR w/ Global Measurements, Complete Interval Values & ST information.



- Capture Unlimited ECG / EKG Durations of Patient ECG Data to Hard Drive, Floppy or any Storage Medium.
- USB or Notebook (PCMCIA) Formats.
- User Selectable Lead Formats on Final ECG Report.
- Print Final Report to Non-fading Laser Hard Copy.
- WINDOWS: Point & Click User Interface.




- Real-Time ECG / EKG Monitoring; 1, 2, 3, 6 or 12 User Selectable Leads
- Complete Test Storage of ECG / EKG Data to Hard Disk, Floppy or Optical.
- 29 Preprogrammed Protocols with 5 User Defined Protocols each with 40 Stages.
- Graphs / Tables: H.R. and BP Trend, Speed / Grade, Representative beats for all Stages.
- View / Print Full Disclosure Measurements for all Stages, ST max. / ST avg..
- On-screen Status of Test: Phase Time, Total Time, Speed, Grade, Target HR, 12 Lead ST Levels & Slopes.
- Report Transmission via Fax or Network.





A free [Database Management System](#) is included with all of our products! This integrates histories for all patients, whether stress, resting, blood pressure, or holter. We are currently developing new tests which we will give you more details on soon.

Sample Screen



AEI Technologies, Inc.
201 Hunters Crossing Blvd.; Suite 10-171
Bastrop, TX 78602
800-793-7751
630-548-3545
Fax: 630-548-3546
Sales@aeitechnologies.com
www.aeitechnologies.com