

S-3A OXYGEN ANALYZER

S-3a I

The Model S-3A Oxygen Analyzer continuously measures the concentration of oxygen.

There are two models:

The S-3A/I Single-Channel Analyzer and the S-3A/II Dual-Channel Analyzer.

Both instruments feature the same high accuracy, rapid response, fine sensitivity and exceptionally long life.



S-3a II



Exclusive Features:

- High 0.01% accuracy at all oxygen concentrations
- Rapid 100 millisecond response ideal for accurate breath-by-breath measurements
- High sensitivity plus exceptional stability capable of detecting changes of 0.001% oxygen at air composition
- Solid ceramic oxide sensing element for long term maintenance-free operation
- Single gas calibration

Exact VO₂, VCO₂ & RER Measurements are too important for speculation, choose the Gas Analyzers researchers have depended on for over 45 years.

**CONTACT US at SALES@AEITECHNOLOGIES.COM or + 1 630 548 3545 for a QUOTE
visit us at www.aeitechnologies.com**